

TRIPLE PACKAGING FOR SPECIMEN SHIPMENT

FOR DIAGNOSTIC AND INFECTIOUS SPECIMENS



PRIMARY LEAKPROOF SPECIMEN CONTAINERS

Labeled with appropriate patient information.



BUBBLE WRAP OR OTHER CUSHIONING MATERIAL

Wrap individual primary specimen containers to prevent contact between specimen containers.



ABSORBENT SHEET BETWEEN PRIMARY SPECIMEN CONTAINER AND BIOHAZARD BAG

Amount of absorbent must be enough to absorb contents.



SECONDARY PACKAGING LEAKPROOF BIOHAZARD BAG



TYVEK ENVELOPE, PRESSURE TESTED TO 95 kPa

ALL OF ABOVE GO INTO OUTER BOX FOR SHIPPING

Use cushioning material to reduce specimen movement within the outer container!